




7. Installing the battery




1. Turn off the meter by pressing the “M” button and hold for 2 seconds before replacing the batteries. Slide down the battery cover on the back of the meter by pushing the cover in the direction of the arrow and put the cover aside.



2. Remove the two used batteries.



3. Insert 2 new batteries with correct + / - direction.
(Battery: AAA 1.5V 2 Alkaline 24A LR03)



4. Put the battery cover back in place by pushing the cover in the opposite direction of the arrow to close the cover into position.

Note:

1. It is suggested that batteries need to be replaced when the battery icon appears empty. If you keep using the meter, the battery icon will start flashing and an error message E_3 will appear on the screen.
2. Note that replacing the battery will not reset the stored testing results.
3. You need to reset the time and date after the batteries are replaced.
4. If there is any unexpected symbol on the display, please follow the above procedure to replace the batteries again.
5. Please dispose of batteries according to your local ordinances.